

Software Version : 6.3.1.0504
 Operator : manager
 Sample Number : 011
 AutoSampler : SER200
 Instrument Name : PE2 PDA
 Instrument Serial # : None
 Delay Time : 0.00 min
 Sampling Rate : 2.2727 pts/s
 Sample Volume : 1.000000 ul
 Sample Amount : 1.0000
 Data Acquisition Time : 19/05/2022 17:18:16

Date : 19/05/2022 17:39:04
 Sample Name : bromoamide
 Study :
 Rack/Vial : 0/21
 Channel : A
 A/D mV Range : 1000
 End Time : 20.00 min

Area Reject : 0.000000
 Dilution Factor : 1.00
 Cycle : 2

Raw Data File : C:\HPLC\Ben\Data\BL-125\230nm_original_solid.raw

Result File : c:\hplc\ben\data\bl-125\230nm_original_solid.rst [Editing in Progress]

Inst Method : C:\HPLC\Ben\Methods\Reverse Phase\50% H2O, 50% MeCN to 40% H2O, 60% MeCN, A230 nm, B220 nm, 2.5 uL injection, 10 min (80 mM default) from C:\HPLC\Ben\Data\BL-125\230nm_original_solid.raw

Proc Method : C:\HPLC\Ben\Methods\Reverse Phase\50% H2O, 50% MeCN to 40% H2O, 60% MeCN, A230 nm, B220 nm, 2.5 uL injection, 10 min (80 mM default) from c:\hplc\ben\data\bl-125\230nm_original_solid.rst [Editing in Progress]

Calib Method : C:\HPLC\Ben\Methods\Reverse Phase\50% H2O, 50% MeCN to 40% H2O, 60% MeCN, A230 nm, B220 nm, 2.5 uL injection, 10 min (80 mM default) from c:\hplc\ben\data\bl-125\230nm_original_solid.rst [Editing in Progress]

Report Format File: C:\HPLC\Ben\Methods\Reverse Phase\50% H2O, 50% MeCN to 40% H2O, 60% MeCN, A230 nm, B220 nm, 2.5 uL injection, 10 min (80 mM default).rpt

Sequence File : C:\HPLC\Ben\Sequences\Sequence Archive\BL-177 Remaining Samples 2.seq

DEFAULT REPORT

Peak #	Time [min]	Area [μV·s]	Height [μV]	Area [%]	Norm. Area [%]	BL	Area/Height [s]
1	3.513	932006.30	121622.93	78.57	78.57	*MM	7.6631
2	3.975	254243.73	30887.47	21.43	21.43	*MM	8.2313
		1186250.03	152510.40	100.00	100.00		

Missing Component Report

Component Expected Retention (Calibration File)

All components were found