



All experiments were carried out on a Thermo Scientific Q Exactive FTMS, employing ASAP / APCI probes.

[M]<sup>+</sup>

m/z	Theo. Mass	Delta (ppm)	Composition
244.1823	<b>244.1822</b>	<b>0.42</b>	<b>C17 H24 O</b>

[M+H]<sup>+</sup>

m/z	Theo. Mass	Delta (ppm)	Composition
245.1900	<b>245.1900</b>	<b>-0.05</b>	<b>C17 H25 O</b>