

EMR6-806-NMR2 - product after column?

Av-400A

Current Data Parameters

NAME Jun22-2017
EXPNO 110
PROCNO 1

F2 - Acquisition Parameters

Date_ 20170622
Time 14.06
INSTRUM spect
PROBHD 5 mm PABBO BB/
PULPROG zg30
TD 32768
SOLVENT CDCl3
NS 16
DS 2
SWH 8012.820 Hz
FIDRES 0.244532 Hz
AQ 2.0447233 sec
RG 126.47
DW 62.400 usec
DE 6.50 usec
TE 296.0 K
D1 1.00000000 sec
TD0 1

===== CHANNEL f1 =====

SFO1 400.4124025 MHz
NUC1 1H
P1 10.90 usec
PLW1 19.95299911 W

F2 - Processing parameters

SI 16384
SF 400.4100000 MHz
WDW EM
SSB 0
LB 0.55 Hz
GB 0
PC 1.00

